

Medical Absorbent Mat

MATE1130 Heavyweight, 90cm x 50m, 1 roll

Prevent slips and falls in the operating theatre with super-absorbent Medical Mat.



- Keeps floors dry — poly-backing keeps floors safe and dry during wet operating procedures
- Ideal for hospital floors — contain and absorb spills from procedures like knee and hip operations where saline fluids are used in abundance for irrigation
- Super absorbent — made of a non-wetting top layer to keep the standing surface dry, a super-absorbent central layer and a non-slip backing that grips to the floor
- Reduce slips and falls — place with white, floor-gripping side down to add traction and reduce slip hazards
- Minimise tripping hazards — no more tripping over bunched-up blankets; increases focus on the patient and not the work environment
- Reduces clean up time — can be easily lifted and disposed of after use
- Simple to use – cut roll to length to fit your space; place mat on clean, dry floor with white, floor-gripping side down to catch liquids as they accumulate; when mat is saturated, lift from floor and remove for disposal in accordance with local regulations
- Heavy-weight mat is only 3.19mm thick

Colour	Medical blue with white base layer
Dimensions	90cm W x 50m L
Mat Weight	Heavyweight
Fluids Absorbed	Water-Based Liquids Only
Absorbency	Up to 166.5 L per bag
Perforated	Not Perforated
Sold as	1 roll per bag
Weight	14.4 kg
# per Pallet	16
Composition	Top layer – Three-dimensional perforated LDPE;

Core - Air-laid cellulose with polypropylene;
Base layer – Polyethylene high-grip film

Packaging	Bag
Poly-Backing	Yes
UV Resistant	No

Technical Documents

Storage & Shelf Life of PIG Absorbents & Absorbent Spill Kits



Hogs Hill • Watt Place • Hamilton International Technology Park • Blantyre • Glasgow • G72 0AH
0800 919 900 • newpig.co.uk • pigpen@newpig.com